





# Mandibular Advancement Splint Product Comparison

	<b>Adjustable PM Positioner</b> 	<b>EMA</b> 	<b>OASYS</b> 	<b>SomnoDent</b> 
Appliance holds jaw securely in place	X			X (must specify)
Allows for speech		X	X	X
Allows for drinking		X	X	X
Allows lips to close	X	X	X	X
Opens nasal passages			X	
Appropriate for a bruxer			X	X (can be customized)
Appropriate for edentulous patients	X (upper only)		X (upper only)	X (upper only)
Adjustable by patient		X		X
CPAP Compatible	X (must specify)		X	X (must specify)
Warranty	90 days		1 year	3 years

## **Contra-indications for mandibular advancement splints**

- Loose teeth
- Advanced periodontal disease
- Central sleep apnea
- Under the age of 18